Instructions for Use and Operation of  
**ColPaC® MASTER CHILLING UNITS**

**C-2 • C-5 • C-6**

ISO 13485/ ISO 9001 CERTIFIED

# TABLE OF CONTENTS

ENGLISH .....	1
DEUTSCH .....	9
FRANÇAIS .....	17
ESPAÑOL .....	23
ITALIANO .....	31
PORTUGUÊS .....	39
NEDERLANDS .....	47
DANSK .....	55
SVENSKA .....	63
PARTS LISTS .....	71

# ENGLISH

## Table of Contents

12 Month Warranty . . . . .	2
Foreword . . . . .	2
Precautionary Instructions . . . . .	3
Instructions for Patient Care and Comfort . . . . .	3
General Information . . . . .	4
General Operating Instructions . . . . .	4
Installation . . . . .	5
Maintenance . . . . .	5
Trouble Shooting . . . . .	6
Technical Specifications . . . . .	7

# 12 Month Warranty

Chattanooga Group ("Company") warrants that the ColPaC® C-2, C-5, and C-6 Master Chilling Unit(s) ("Product") is free of defects in material and workmanship. This warranty shall remain in effect for one (1) year from the date of original consumer purchase of this Product and extends to any owner of the Product during the warranty period. If this Product fails to function during the one year warranty period because of a defect in material and/or workmanship, Company or the selling dealer will replace or repair this Product without charge within a period of 30 days from the date on which the defective Product is returned to the Company or the dealer. Company or the dealer will ship the replacement or the repaired Product to the consumer's residence.

## This Warranty Does Not Cover

1. Replacement parts or labor furnished by anyone other than the Company, the dealer or an approved Company service agent.
2. Defects or damage caused by labor furnished by someone other than Company, the dealer or an approved Company service agent.
3. Any malfunction or failure in the Product while it is in the possession of the owner during the warranty period if the malfunction or failure is not caused by a defect in material or workmanship, or if the malfunction or failure is caused by unreasonable use, including the failure to provide reasonable and necessary maintenance.

## Company Shall Not Be Liable for Incidental or Consequential Damages to Property or Business

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

TO OBTAIN SERVICE from Company or the selling dealer under this warranty, the owner must do or abide by the following:

1. A written claim must be made within the warranty period to the Company or the selling dealer. If the claim is made to the Company, written claim should be sent to:

Chattanooga Group  
4717 Adams Road  
P.O. Box 489  
Hixson, TN 37343 USA  
Phone: +1 (423) 870-7200 Outside U.S.A.  
Fax: +1 (423) 870-2046

2. The Product must be returned to the Company or the selling dealer by the owner.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state.

The Company does not authorize any person or representative to create for it any other obligation or liability in connection with the sale of the Product. Any representation or agreement not contained in the warranty shall be void and of no effect.

## Foreword

This manual has been prepared for the owners and operators of the ColPaC Models C-2, C-5, and C-6 Master Chilling Units. It contains instruction on operation, precautions, maintenance and parts information. To obtain maximum life and efficiency from your ColPaC C-2, C-5, and C-6 units, and to aid in its proper operation, read and understand this manual thoroughly and follow all instructions prior to operating the unit.

The specifications put forth in this manual were in effect at the time of publication. However, owing to Chattanooga Group's policy of continuous improvement, changes to these specifications may be made at any time without obligation on the part of Chattanooga Group.

## Precautionary Instructions

1. **CAUTION:** Read, understand and practice the Precautionary and Operating Instructions. Know the limitations associated with the ColPaC Master Chilling Units. Observe the cautionary and operational decals installed on the units.
2. **CAUTION:** Do not place chemical activated and/or instant freeze type packs in the ColPaC Master Chilling Unit. This type of pack may break and cause corrosive damage to the freezer compartment and refrigeration system.
3. **WARNING:** Make certain that the unit is electrically grounded by connecting to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.

## Instructions for Patient Care and Comfort

1. Know the cooling characteristics of the ColPaC.
2. Exercise extreme care when applying the ColPaC.
3. Know your patient: Use extreme care when applying ColPaCs to areas where the patient may have loss of feeling.
4. Always instruct the patient to inform you if the ColPaC is too cold.
5. Always place one or more layers of dampened, thin towel between the skin and the ColPaC if the application becomes painful. A warm, damp towel reduces the initial shock of the cold application.
6. Use care and judgment when applying ColPaCs to infants and elderly. Constantly check the temperature of the ColPaCs.
7. Exercise extreme caution when applying ColPaCs to patients with thin skin or to "bony" areas: shoulders, elbows, knees, etc.
8. Do not use the ColPaC in the presence of poor circulation or open wounds.
9. Carefully monitor the first treatment with the ColPaC and adapt the method of application to the tolerance and needs of the patient.
10. Tolerance of cold, like tolerance of heat, will vary from patient to patient.
11. During treatment, frequently check the condition of the patient and the area being treated.
12. Patient should not sit, lie or sleep on the ColPaC.
13. Treatment time should not exceed 30 minutes.

## General Information

Listen carefully after connecting the unit to the electrical power to make sure the compressor goes into operation immediately. A light frost covering should form on the dividers in the chilling chamber. Insert the individual ColPaCs between the dividers in the chilling chamber and close the lid. The various ColPaCs are shown below.

The initial ColPaC chilling, starting from room temperature, will require several hours. The time needed to re-chill the ColPaC after use on a patient will depend on the length of time in use with the patient and the time out of the unit. Normally a short time is required for re-chilling your ColPaCs.

**Note:** While in operation, the top of the unit will become warm because of a specially designed heating system installed under the top rim. This is to eliminate condensation of moisture around the top of the unit.

### Care of ColPaCs

Always return the ColPaC to the chilling chamber after each treatment where it is kept chilled, clean and ready for immediate use. The ColPaC may be cleaned by scrubbing it with soap and water. It is best not to fold the ColPaC for storage. If the outside plastic cover of the ColPaC should become torn or cut, minor repairs may be made with plastic or vinyl tape. Clean the surface carefully and apply tape at room temperature. Always replace ColPaCs when the plastic covering becomes damaged beyond repair. Original ColPaCs may be ordered from an authorized Chattanooga Group dealer.

The temperature of the ColPaC, when chilled, is approximately 10° F to 21° F (-6° C to -12° C). Be familiar with the cooling characteristics of the ColPaC.

## General Operating Instructions

The ColPaC Master Chilling Unit is simple and easy to operate. The thermostat controls the temperature of the chilling chamber. It is located in the rear of the unit, behind the removable panel. The refrigeration system of the ColPaC Master Chilling Unit has been adjusted, tested and calibrated to chill the ColPaCs to a temperature between 10° F to 21° F (-6° C to -12° C). Extensive clinical testing and observation have shown this to be the optimum temperature range. It is recommended that this adjustment not be altered until the ColPaC unit and the ColPacs have been used for some time, and then only if adjustment of the chilling temperature is absolutely necessary. A slight adjustment can be made by removing the rear panel and turning the adjustment knob on the thermostat. To lower the temperature (colder), the knob should be turned clockwise. To raise the temperature (warmer), the knob should be turned slightly counter-clockwise.

**CAUTION:** The thermostat is extremely sensitive and the slightest adjustment will alter the temperature several degrees. The temperature of the chilling chamber should be checked with a thermometer after every adjustment, before using the ColPaCs.

The ColPaC Master Chilling Unit operates through the use of a hydraulic capillary type thermostat and can be left connected to the power source, providing a ready supply of chilled ColPaCs. Locate the unit away from humid areas to avoid excessive frost build-up on the refrigeration coils.

## General Operating Rules

- Always keep the ColPaCs in the chilling chamber of **the unit when not in use**. This keeps the packs clean, chilled, and ready for use, and prevents pack damage.
- Never place a damaged ColPaC into the chilling chamber of the unit.
- Check for ice/frost build-up daily and defrost unit when ColPaCs become crowded in chamber.
- Never adjust the thermostat too high or too low. See caution above.
- Always replace ColPaCs when the plastic covering becomes damaged beyond repair or when ColPaC becomes firm or hard and not moldable around the area being treated. Minor repairs may be made with plastic or vinyl tape.
- If the unit is to be left unattended for an extended period of time, unplug and follow the defrosting procedures found in the "Maintenance" section of this manual
- Always use factory authorized replacement parts.
- Always clean and sterilize ColPaCs, which are exposed to contamination. Do not sterilize ColPaCs with boiling water. Gas sterilizing is not recommended. Alcohol or similar substance may be used to decontaminate the packs. If any contamination is suspected, destroy the ColPaCs.

## Installation

Remove all packing material inside the ColPaC chilling chamber. It is important to remove all loose pieces of paper and other packing material, which might otherwise plug the drain! Make certain that the drain valve is closed.

Check the voltage rating on the serial plate which is located on the back of the unit. Plug the unit into properly grounded outlet/power source. Do not attempt to use Direct Current (D.C.). Follow the procedures indicated in the "Precautionary Instructions." Do not attempt to use the unit if it is not properly grounded.

### Explanation of Symbols on Equipment

Attention: See accompanying documentation

Type B Equipment: An adequate degree of protection against electric shock is provided, particularly regarding allowable leakage currents and reliability of protective earth connection.

## Maintenance

### Care of Unit

The ColPaC Master Chilling Unit is equipped with a hydraulic capillary type thermostat in order to maintain the proper ColPaC temperature evenly in the chilling chamber. The unit may be left connected to the electrical outlet, providing a ready supply of chilled ColPaCs.

Periodically, the accumulated frost in the chilling chamber will build up to a point where the ColPaCs are crowded and insertion or removal of the packs becomes difficult. The units should be defrosted when this occurs and the inside of the tank wiped dry with a clean towel or cloth.

### Defrosting Unit

1. Disconnect the unit from its power source.
2. Raise the chamber lid and remove ColPaCs. Leave lid open to speed defrosting. A small amount of warm water may be poured into the unit to speed defrosting.
3. Attach extension drain hose to the valve, located beneath chilling chamber, behind rear access panel.
4. Open bottom drain valve.
5. Allow unit to defrost completely. All excess water must be removed from the chilling chamber before the unit is started again.

6. Close bottom drain valve and detach extension drain hose.
7. Plug unit back into the power source.
8. Clean and return ColPaCs to the chilling chamber.
9. Close chamber lid - chilling sequence begins again.

Dust and lint will accumulate on the condenser fins of the refrigeration unit. This must be removed periodically to ensure highest efficiency. The condenser may be reached by removing the rear access panel. Always disconnect power before cleaning the condenser.

### Replacing Thermostat

1. Remove plate.
2. Remove bushing nut.
3. Remove packing in bushing nut.
4. Remove brass nut.
5. Remove rubber washer.
6. To reinstall, reverse procedure.
7. After installation, seal top of packing nut with silicone sealant.

## Troubleshooting

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTION
A. Unit plugged in but packs not chilled.	<ol style="list-style-type: none"> <li>1. Breaker for electric circuit off</li> <li>2. Thermostat not functioning</li> <li>3. Heating element burned out</li> <li>4. Compressor failure</li> <li>5. Refrigerant leak.</li> </ol>	<ol style="list-style-type: none"> <li>1. Turn breaker ON</li> <li>2. Replace thermostat</li> <li>3. Replace element</li> <li>4. Contact factory</li> <li>5. Contact factory.</li> </ol>
B. Packs too cool.	<ol style="list-style-type: none"> <li>1. Thermostat set too low.</li> <li>2. Thermostat failure.</li> </ol>	<ol style="list-style-type: none"> <li>1. Adjust thermostat.</li> <li>2. Replace thermostat.</li> </ol>
C. Packs too warm	<ol style="list-style-type: none"> <li>1. Thermostat set too high</li> <li>2. Thermostat failure</li> <li>3. Power failure.</li> </ol>	<ol style="list-style-type: none"> <li>1. Adjust thermostat</li> <li>2. Replace thermostat</li> <li>3. Check electrical circuit.</li> </ol>

## Technical Specifications (115 Volt)

<b>Rated Voltage:</b>	115V
<b>Rated Frequency:</b>	60 Hz
<b>Current:</b>	4 Amps
<b>Work:</b>	.14914KWh or 1/5 hp
<b>Type of Refrigerant:</b>	R-134A
<b>Amount of Refrigerant:</b>	C-2, 198 g (7 oz.) C-5, 198 g (7 oz.) C-6, 213 g (7.5 oz.)
<b>Test Pressures, All Models:</b>	High, 1.62 MPa (235 psi) Low, 1.03 MPa (150 psi)
<b>Normal Operating Temperature,</b>	
<b>    All Models:</b>	-12° C to -6° C (+10° F to +21° F)

## Technical Specifications (220 Volt)

<b>Rated Voltage:</b>	220-240V~
<b>Rated Frequency:</b>	50 Hz
<b>Current:</b>	2 Amps
<b>Circuit Breaker:</b>	10 Amp, Anti-Surge
<b>Work:</b>	.14914KWh or 1/5 hp
<b>Type of Refrigerant:</b>	R-134A
<b>Amount of Refrigerant:</b>	C-2, 198 g (7 oz.) C-5, 198 g (7 oz.) C-6, 213 g (7.5 oz.)
<b>Test Pressures, All Models:</b>	High, 1.62MPa (235 psi) Low, 1.03 Mpa (150 psi)
<b>Normal Operating Temperature,</b>	
<b>    All Models:</b>	-12° C to -6° C (+10° F to +21° F)

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I

## Parts Lists Model C-2 (220 volt)

Item	Part No.	Description	Qty
1	20500	THERMOSTAT, Resale	1
2	22146	LAMP, Indicator	1
3	21283	CORD SET, w/J32	1
4	21635	HOSE, Nylon, 5/8" x 12, Coupled, M-F	1
5	21983	BACK DOOR, Sub-Assembly	1
6	20514	HANDLE, Lid, Resale	1
7	22219	LID, Assembly	1
8	21373	VALVE, Boiler Drain, ½" ID, Chrome	2
9	21259	CASTER, 3" EY4520, 7/16-20 x ½"	4
10	22122	COMPRESSOR, 1/5 HP, BM2A529-1	2
11	4228	THERMOMETER, Dial	2
12	22409	COIL COVER, 13.5" X 22.75", Outside	1
13	22420	COIL COVER, Sub-Assembly	1
14	22955	COIL COVER, Middle Aluminum	3
15	20237	HANDLE, Decorative, 12"	2
16	21011	SWITCH, DPST, 2GK51-NB23, Common	1
17	22173	BREAKER, Heinemann, KD1-10	2
18	60286	LATCH, Door, Rear	1

## Parts List Model C-2 (115 volt)

Item	Part No.	Description	Qty
1	20500	THERMOSTAT, Resale	1
2	22148	LAMP, Pilot Light	1
3	21282	CORD SET, w/J28G	1
4	21635	HOSE, Nylon, 5/8" x 12", Coupled, M-F	1
5	21983	BACK DOOR, Sub-Assembly	1
6	21260	CASTER, 3" G421-3	1
7	20514	HANDLE, Lid, Resale	1
8	22219	LID, Assembly	1
9	21373	VALVE, Boiler Drain, ½" ID, Chrome	1
10	21259	CASTER, 3" EY4520, 7/16-20 x ½"	4
11	22121	COMPRESSOR, BM2A528-1	1
12	4228	THERMOMETER, Dial	1
13	22409	COIL COVER, 13.5" X 22.75", Outside	1
14	22420	COIL COVER, Sub-Assembly	1
15	22955	COIL COVER, Middle, Aluminum	3
16	20237	HANDLE, Decorative, 12"	2
17	60286	LATCH, Door, Rear (not shown)	1

## Parts Number Model C-5 (220 volt)

Item	Part No.	Description	Qty
1	20500	THERMOSTAT, Resale	1
2	22146	LAMP, Indicator	1
3	21283	CORD SET, w/J32	1
4	21635	HOSE, Nylon, 5/8" x 12, Coupled, M-F	1
5	22457	BACK DOOR, Sub-Assembly	1
6	20514	HANDLE, Lid, Resale	1
7	22297	LID, Assembly	1
8	21373	VALVE, Boiler Drain, ½" ID, Chrome	1
9	21259	CASTER, 3" EY4520, 7/16-20 x ½"	4
10	22122	COMPRESSOR, 1/5 HP, BM2A529	1
11	4228	THERMOMETER, Dial	1
12	22461	COIL COVER, Outside	1
13	22462	COIL COVER, Sub-Assembly	1
14	22956	COIL COVER, Middle Aluminum	3
15	20237	HANDLE, Decorative, 12"	2
16	21011	SWITCH, DPST, 2GK51-NB23, Common	1
17	22173	BREAKER, Heinemann, KD1-10	2
18	60286	LATCH, Door, Rear	1

## Parts List Model C-5 (115 volt)

Item	Part No.	Description	Qty
1	20500	THERMOSTAT, Resale	1
2	20042	LAMP, Pilot Light	1
3	21282	CORD SET, w/J28G	1
4	21635	HOSE, Nylon, 5/8" x 12", Coupled, M-F	1
5	22457	DOOR, Rear, Complete	1
6	21260	CASTER, 3" G421-3	1
7	20514	HANDLE, Lid, Resale	1
8	22297	LID, Assembly	1
9	21373	VALVE, Boiler Drain, ½" ID, Chrome	1
10	21259	CASTER, 3" EY4520, 7/16-20 x ½"	4
11	22121	COMPRESSOR, BM2A528-1	1
12	4228	THERMOMETER, Dial	1
13	22461	COIL COVER, Outside	1
14	22462	COIL COVER, Sub-Assembly	1
15	22956	COIL COVER, Middle, Aluminum	3
16	20237	HANDLE, Decorative, 12"	2
17	60286	LATCH, Door, Rear (not shown)	1

## Parts Number Model C-6 (220 volt)

Item	Part No.	Description	Qty
1	20500	THERMOSTAT, Resale	1
2	22146	LAMP, Indicator	1
3	21283	CORD SET, w/J32	1
4	21635	HOSE, Nylon, 5/8" x 12, Coupled, M-F	1
5	22457	BACK DOOR, Sub-Assembly	1
6	20514	HANDLE, Lid, Resale	1
7	22297	LID, Assembly	1
8	21373	VALVE, Boiler Drain, ½" ID, Chrome	1
9	21259	CASTER, 3" EY4520, 7/16-20 x ½"	4
10	22122	COMPRESSOR, 1/5 HP, BM2A529	1
11	4228	THERMOMETER, Dial	1
12	22461	COIL COVER, Outside	1
13	22462	COIL COVER, Sub-Assembly	1
14	22956	COIL COVER, Middle Aluminum	3
15	20237	HANDLE, Decorative, 12"	2
16	21011	SWITCH, DPST, 2GK51-NB23, Common	1
17	22173	BREAKER, Heinemann, KD1-10	2
18	60286	LATCH, Door, Rear	1

## Parts List Model C-6 (115 volt)

Item	Part No.	Description	Qty
1	20500	THERMOSTAT, Resale	1
2	20042	LAMP, Pilot Light	1
3	21282	CORD SET, w/J28G	1
4	21635	HOSE, Nylon, 5/8" x 12", Coupled, M-F	1
5	22457	DOOR, Rear, Complete	1
6	21260	CASTER, 3" G421-3	1
7	20514	HANDLE, Lid, Resale	1
8	22297	LID, Assembly	1
9	21373	VALVE, Boiler Drain, ½" ID, Chrome	1
10	21259	CASTER, 3" EY4520, 7/16-20 x ½"	4
11	22121	COMPRESSOR, BM2A528-1	1
12	4228	THERMOMETER, Dial	1
13	22461	COIL COVER, Outside	1
14	22462	COIL COVER, Sub-Assembly	1
15	22956	COIL COVER, Middle, Aluminum	3
16	20237	HANDLE, Decorative, 12"	2
17	60286	LATCH, Door, Rear (not shown)	1

23620B

© 2002 Encore Medical

ISO 13485/ ISO 9001 CERTIFIED

4717 Adams Road  
P.O. Box 489  
Hixson, TN 37343 USA  
Phone: (423) 870-7200  
Fax: (423) 870-2046  
[www.chattgroup.com](http://www.chattgroup.com)